



Addendum EUROLineScan






Appendix	1
EUROLineScan	1
1 ImpEx.....	3
1.1 General information	3
1.2 Language settings	3
1.3 Export & Import.....	3
1.3.1 Worksheets/ Exporting lots.....	3
2 Online connection	7
2.1 GDT Connection.....	7
2.1.1 Settings.....	7
2.1.2 Importing protocol configurations/runs	7
2.1.3 Exporting the evaluation results	8
2.1.4 Tests.....	10
2.2 ELO connection (EUROLabOffice).....	11
2.3 Word connection.....	11
3 Client Server Installation	12
4 Quick Data Source Changeover	14



1 ImpEx

1.1 General information



The ImpEx program is provided together with EUROLineScan program and can normally be found in the subdirectory „ImpEx“ of the EUROLineScan program. The ImpEx program was designed to enable the exporting of lot or protocol data from another database of the EUROLineScan program and to later reimport. The protocol data and lots are XML-Format secured and are automatically packed into a zip. file. When importing, a zip. archive can be selected enabling the included data in that program to then be imported into another selected database.

1.2 Language Settings


The ImpEx program can be executed in German and English. To change the language select **Einstellungen → Sprache einstellen** or **Settings → Set language** and select the desired language. Restart the program to activate the selected language.

1.3 Export & Import

Before data is exported or imported, please check if the correct database, user name and password has been set.

1.3.1 Protocols/ Exporting lots


To export protocols or lots from the EUROLineScan database, please click on the corresponding button **Protokolle exportieren/export protocols** or **Chargen exportieren/export lots**. Enter the file name and the directory where the data is to be stored. Usually the export data is stored in the EUROLineScan installation directory: „ImpEx\Archiv“.



- To store the exported data on your *Desktop*, select *Desktop* first then enter the desired file name.
- To store the exported data on a removable disk (**USB-Stick**), first select *Arbeitsplatz /My computer* then select removable disk and enter the file name.

Now click to save. In the following window you can select the **protocols and/or lots** that are to be exported. Input fields can be filled on the right side which enables only selected protocols and lots to appear in the list. Please give attention not only to the protocols that are to be exported, but also the the **export settings**:

- **Export existing protocol images**
When choosing this setting, the protocol images that were scanned in and saved in the EUROLineScan program will be exported.
- **Export patient data**
When choosing this setting, data such as patient name, birth date, lab etc. will be exported with the protocol. Should the patient data NOT be exported, then deactivate this setting.






Click Ok to start the process. During this process you will see the following window. This process may take a while depending on how many protocols or lots have been selected.



As soon as the process is completed you will be informed if the export was successful or not. Should an error occur an export file will not be created.

1.3.2 Importing protocols and lots


To import protocols and lots from an export file back into the EUROLinScan database click on the respective button **Import protocols** or **Import lots**. In the next window enter in the export file of which data is to be imported. Usually the export data is stored in the EUROLinScan installation directory: *VmpEx\Archiv*, for this reason the program initially offers the archive directory for selection of the export file.




- To open an export file stored on your **Desktop**, first select *Desktop* and then the respective file.
- To open an export file stored on an removable disk (**USB-Stick**), first select *Arbeitsplatz/My computer*, then removable disk and the respective file.



Click on open to start importing the data. The following window will appear during the import. This process may take a while depending on how many protocols or lots have been selected.



As soon as the process is completed you will be informed if the export was successful or not. Should an error occur the import up to that point will be cancelled, so that the EUROLinScan database is in the same state as it was in before the import process began. A list of all the imported data is shown if the import was successful.





2 Online connection

The EUROLinScan program can exchange data with other programs, concerning the compilation of a protocol (import) and the issuing of evaluation results (export).

On the „Connection“ page you can select which option is to be used. On the sub-pages further settings can be made to each option.

2.1 GDT connection

Data exchange occurs through individual files for each strip of the protocol. The data is based on a GDT (Gerätedaten-Träger) format. A detailed interface description is available on the website: <http://www.qms-de.org/gdt/gdtenglisch.zip>.

2.1.1 Settings




Illustration 1 GDT settings in EUROLinScan

The settings for the GDT connection can be entered in the directory for the import („Source directory“) and export („Destination directory“). These directories cannot be the same. If „test ident without test name“ has been chosen, the test will be left out before the band name in field 8410.

To skip the label, alignment bar and indicator bands when listing the bands the checkbox „band results without label and indicator“ must be selected. Furthermore, when selecting „do not check the length“ a warning can be disabled during verification of the field 8100.

For allergy strips the concentration levels can be entered instead of the intensity when „if possible use concentration instead of intensities“ is selected.

2.1.2 Importing protocol configurations

All data in the import directory with the ending „gdt“ are sorted and read alphabetically. After reading, the files are deleted by the EUROLinScan program. A file must exist for each strip of the protocol. Record type 6302 is used („New Test Request“).

eg. „unter000.gdt“ (field identification in bold type):

```
01380006302
014810000116
014921802.10
014300012345
0193101Mustermann
0143102Peter
017310313012000
0118410BG
```




The fields 8000, 8100, 9218 are standard fields and are sufficiently described in the GDT format documentation.

The EUROLinScan program requires the patient ID for field 3000. This will clearly identify the patient in the database. This field can contain up to 64 characters as opposed to the official GDT format.

Also in the fields 3101 and 3102 for the last name and first name of the patient 64 characters are available. Patient birth date is entered into field 3103.

A test abbreviation is required from the test list for field 8410. (See chapter 2.1.4)

Furthermore the following fields can be entered:

Field	Content
6330	Contains the character string „Labornummer“
6331	Alphanumerical content for the lab number with max. 32 characters
6332	Contains the character string „Einsender“
6333	Alphanumerical content for “Sent from” with max. 32 characters
6334 oder 8300	Contains the character string “Einsenderdatum”
6335 oder 8301	Date of receipt sample in format DDMMYYYY
6336	Contains the character string „SID“
6337	Alphanumerical content for test kit with max. 32 characters

2.1.3 Exporting the evaluation results

For each strip of the protocol a file is created with the file format name „unter000.gdt“. The numbers „000“ represent successive numbering that are chosen to prevent overwriting of the existing files. Deleting files must be done manually or with a different software.

The evaluation for EUROLIN and Westernblot must be differentiated. For the EUROLIN each individual band represents one result, whereas in Westernblot only from all of the bands can one or several results be concluded.

Record type 6310 is used („Test Data Transfer“).

eg. „unter000.gdt“ (field identification in bold type):

```

01380006310
014810000936
014921802.10
0210103EUROLinScan
022300021009Ea 30µl
0158402ALLG00
017620012082004
0156201125417
0166330PRT/STR
01563311534/0
0118410VG
0338411Borrelia EUROLIN-WB IgG
017843212082004
0158439125417
0168480positive
0128410Al
0218411Alignment bar
01084200
0198421Intensität
0108480o
0128410p17
0138411p 17

```




011842024
0198421Intensität
0108480+
0128410p19
0138411p 19
011842034
0198421Intensität
0108480+
0128410p21
0138411p 21
0128420143
0198421Intensität
0108480+
0128410p25
0208411p 25, Osp C
01084201
0198421Intensität
0108480o
0128410p30
0138411p 30
0128420128
0198421Intensität
0108480+
0128410p31
0208411p 31, Osp A
011842051
0198421Intensität
0108480+
0128410p39
0208411p 39, Bmp A
0128420159
0198421Intensität
0108480+
0128410p83
0138411p 83
011842084
0198421Intensität
0108480+
0138410VIsE
0138411VIsE
0128420111
0198421Intensität
0108480+

The fields 8000, 8100, 9218 are standard fields and are sufficiently described in the GDT format documentation.

In field 0103 „EUROLineScan“ is assigned as the software. The PatientID with the field identification 3000 then follows.

„ALLG00“ is assigned for equipment and procedure specification identity (field identification 8402).

For the day and time required for data collection (6200 & 6201) the establishment date of the protocol is used.

Vacant categories follow thereafter. The combination 6330 and 6331 define the internal program IDs for the protocol and the strip number (format: protocol number/strip number). All evaluated results of one protocol must have the same protocol number. The count of the strip numbers of each protocol begins with 0 and increases consecutively. The second vacant category (6332 and 6333) is used for the lab number.

Possible existing results of the strips will follow next. (field 8480).

For the Test-Ident a test abbreviation is required from the test list for field 8410. (See chapter 2.1.4)



In field 8411 the description taken from the list is entered. The incubation date of the protocol is used for the acceptance date and time (8432 & 8439). One or several fields with results will follow.

Thereafter a list with all of the bands of the strip will follow.

In the fields 8410 and 8411 the short and long syntax of the band is included. The measured intensity for these bands is given as a whole number from 0 to 255 in field 8420. The unit used is "Intensität" (field identification 8421). For the result text in field 8480 the symbol from the details view is used.

2.1.4 Tests

Abbreviation	Description
AP3a	Ana Profile3 (13 Ag, SS-A and Ro52)
APa	Ana Profile 1 (10 Ag, SS-A and Ro52)
At	Allergy EUROLINE Atopy
Cr	Allergy EUROLINE Cross reactions
E1a	Ena ProfilePlus 1 (6 Ag, SS-A and Ro52)
EBV	EBV Profile 1 EUROLINE
EBV2	EBV Profile 2 (5 Antigens) EUROLINE
FO	Allergy EUROLINE Food
In	Allergy EUROLINE Inhalation
Ins	Allergy EUROLINE Insects
Pd	Allergy EUROLINE Pediatric Inhalation
TOR	T.O.R.C.H. Profiel
BG	Borrelia burgdorferi IgG
BM	Borrelia burgdorferi IgM
YA	Yersinia enterocolitica IgA
YG	Yersinia enterocolitica IgG
PA	Helicobacter pylori IgA
PG	Helicobacter pylori IgG
GG	Borrelia garinii IgG
GM	Borrelia garinii IgM
HG	Herpes simplex Virus IgG
HM	Herpes simplex Virus IgM
EG	EBV Westernblot IgG
EM	EBV Westernblot IgM
AG	Treponema pallidum IgG
AM	Treponema pallidum IgM
FG	Borrelia afzelii IgG
FM	Borrelia afzelii IgM
KG	Echinococcus granulosus IgG
KM	Echinococcus granulosus IgM
CG	Cytomegalovirus IgG
CM	Cytomegalovirus IgM
VG	Borrelia EUROLINE-WB IgG
VM	Borrelia EUROLINE-WB IgM
Leb	Leber Profile
Med1	Allergy EUROLINE Mediterranean Inhalation
AP3b	Ana Profile3 (14 Ag, PCNA)
APb	Ana Profile 1 with nucleosomes
Blnk	a blank Westernblot strip
ANCA	MPO_PR3 EUROLINE
FO2	Allergy EUROLINE Korea Food2
ISEA	Allergy EUROLINE Inhalation SEA
FSEA	Allergy EUROLINE Food SEA
Ped	Allergy EUROLINE Pediatric
Med3	Allergy EUROLINE Mediterranean Inhalation
IN1	Allergy EUROLINE Korea Inhalation



2.2 ELO connection (EUROLabOffice)

When using the EUROLabOffice program from EUROIMMUN data can be exchanged in a special XML format. If „EUROLabOffice“ is set for the connection an import file will be requested when establishing a new protocol. All information will be taken out of this file to create the protocol (except for the lot for Westernblot tests). The file selection can also be interrupted and the protocol can be entered manually. In the setting for the ELO connection a path can be pre-set with the file name („import file“).

To export („export path“) a path is defined to later direct the file containing the evaluation results. The file name is automatically selected „export[number].xml“ so that no other files are overwritten.

2.3 Word connection

If the setting „Word“ has been selected for the connection, the EUROLinScan program will try to open the designated word file for the word connection settings through „export“. This word file can be freely formulated by the user. For information that is to be inserted from the EUROLinScan text markers must be set at these positions according to the following scheme.

Name of the text mark	Content
strnr	Strip number
unterid	Test ID from the database
unterkuerzel	Test abbreviation
unterkurz	Test short description
unterbez	Test description
patientid	Patient ID
patname	Patient name
patdatum	Patient birth date
labor	Lab number
streifen	Strip number
testkit	Test kit
blutdatum	Blood sample receipt in Lab
bluteinsender	Blood sender
chargeid	Lot ID from the database
charge	Lot description
strerg	Strip results
bnderg	Band results

For each strip of the protocol a word file is opened and filled. At the end the word program will open and the individual files will be displayed.



3 Client server installation

The EUROLinScan program can be installed onto several computers. The program is independently installed and does not exchange data.



Illustration 2 Each set up administrates its own data

If the computers are connected to a network they can use a community data source enabling work and viewing of the protocols and results on different computers.

The data that can be exchanged includes the content of the database and the images of the strips. Additionally the images to be saved for the protocol and templates can be used collectively.




Illustration 3 One computer takes over the administration for all of other computers

This can be easily achieved when one computer is designated as the server and all other set ups can use the data on this computer. This enables a unified data base for all installations.

For this the strips, protocols and templates must be shared to the network. In addition access is not to be blocked from outside on port 3050 by firewalls.

The database itself cannot be released to the network. The interbase database file must always be on the local computer. The EUROLinScan can be installed normally onto the computer which is to serve as a client. A few features can be deselected when installing the interbase database (a part of the EUROLinScan installation) to save resources on the computer.




Illustration 4 Necessary client installation settings



After installation a few changes are necessary. Open the program „BDEADMIN“ in the directory“ C:\Program Files\Common Files\Borland Shared\BDE\“. On the left side select the entry „euroscan“. After this step, on the right side it is possible to edit the entry for „SERVER NAME“. There the IP address of the server computer must be entered. A colon and the path to database file on the server must be entered following this. The standard login for the database is „euroscan“ for the name and password.




Illustration 5 Program BDE administration

In addition, the settings in the EUROLinScan program for the entries „Picture of strips“ , "Protocol pictures“ and „Template pictures“ must be adjusted to the share from the server.




Illustration 6 The administration of data is designated to one server

The third possibility is to save the images and the database to a dedicated server. Another SQL database can be used instead of the interbase database. The program has been successfully tested with Oracle and MySQL databases.

To save images it is not necessary to have a Windows server. Samba under Linux are also successfully in use.



4 Quick Data Source Changeover

In most cases the program is run with one data source. This can be adjusted to the settings on the page "program" with the entries to „Database alias“, „DB-Name“, „DB-Pass“, „Picture of strips“, „Protocol pictures“ and „Template picture“. If an automatic switch between the different data sources is desired, this can be made possible with a special XML file. This XML file can be entered into the program settings „Databases“. The content could look like this example:

```
<databases override="false">
  <database>
    <name>euroscan</name>
    <alias>euroscan</alias>
    <dbname>euroscan</dbname>
    <dbpass>euroscan</dbpass>
    <strips>C:\Programme\EUROIMMUN\EUROLinScan\pictures\strips</strips>
    <templates>C:\Programme\EUROIMMUN\EUROLinScan\pictures\templates</templates>
    <worksheets>C:\Programme\EUROIMMUN\EUROLinScan\pictures\worksheets</worksheets>
    <valid_from>2006-02-02T23:59:59</valid_from>
  </database>
  <database>
    <name>euro2</name>
    <alias>euro2</alias>
    <dbname>euroscan</dbname>
    <dbpass>euroscan</dbpass>
    <strips>\\terra\euroscan\pictures\strips</strips>
    <templates>\\terra\euroscan\pictures\templates</templates>
    <worksheets>\\terra\euroscan\pictures\worksheets</worksheets>
    <valid_from>2005-12-31T23:59:59</valid_from>
  </database>
</databases>
```

In this example two entries in the menu in „database“ would appear after the program has been started. With these two options it is possible to switch between two data sources without having to close the whole program. The individual data sources („database“ elements) are in a „databases“ element. Each „database“ element represents an entry in the program menu. The parameter „override“ shows if the most current data source is to be automatically selected when the program is started.

The „name“ element will take on the description for the data source that is shown in the menu. The elements „Database alias“, „dbname“ and „dbpass“ refer to database and correspond to the entries „database alias“, „DB-Name“ and „DB-Pass“ in the program settings. The directories in „strips“, „templates“ and „worksheets“ correspond to „Picture of strips“, „Protocol pictures“ and „Template pictures“.

With the entry to „valid_from“ in format „YYYY-MM-DDTHH:NN:SS“ (Y=year, M=month, D=day, H=hour, N=minute, S=second) will on the one hand indicate when the entry is valid and on the other also show how recent the entry is in relation to the parameter „override“ (see above).



