

Communication Connections Between the ADVIA® 1650 Workstation and a Host Computer ADVIA® 1650 Chemistry System
Communication Connections
Between the ADVIA® 1650 Workstation
and a Host Computer

Publication No. TN9-9155-16 / December 2000 Supersedes TN9-9155-15, issued November 1999

Publication Change Record

| Publication Number | <u>Date</u> |
|---------------------------|---------------|
| TN9-8155-01 | December 1998 |
| TN9-9155-15 | November 1999 |
| TN9-9155-16 | December 2000 |

List of Effective Pages

| Page Number | <u>Date</u> |
|-------------|---------------|
| Cover | November 1999 |
| Title | December 2000 |
| A-page | December 2000 |
| i | December 2000 |
| 1 | November 1999 |
| 2-3 | December 2000 |
| 4 | November 1999 |
| 5-6 | December 2000 |
| 7 | December 1998 |
| 8-9 | December 2000 |
| 10-11 | December 1998 |
| 12-14 | December 2000 |
| 15 | December 1998 |
| 16 | December 2000 |
| 17 | December 1998 |
| 18 | December 2000 |
| 19 | December 1998 |
| 20-22 | December 2000 |
| 23 | December 1998 |
| 24 | December 2000 |
| 25 | December 1998 |
| 26-30 | December 2000 |
| 31-33 | November 1999 |
| 34 | December 2000 |
| 35 | November 1999 |
| 36-37 | December 1998 |

Contents

| ■ Purpose | 1 |
|---|----|
| ■ Interconnecting the ADVIA 1650 Chemistry System and the Host Computer | 1 |
| Communication Specifications and Parameters | 2 |
| Connector Pin Assignments | 3 |
| Message Frame Format | 4 |
| Communication Modes | 5 |
| Hardware Configuration | |
| System Customization | |
| Available Operating Modes | |
| Sample Identification | 8 |
| Registration Number | 8 |
| Position Number | 8 |
| Sample Number | 8 |
| Sample Information | 8 |
| Classification of ADVIA 1650 Samples | 9 |
| Item Order Operation | 10 |
| Host Item Code Table | 11 |
| Communication Protocols | |
| Example of Response Time to Meet 3.0 Second Limitation | 13 |
| Batch Item Registration Protocol | 14 |
| Automatic Item Registration Protocol | 16 |
| Real Item Registration Protocol | 18 |
| Batch Data Output Protocol | 20 |
| Real Data Output Protocol | 21 |
| Error Processing | |
| Abnormal Text Format | |
| Receipt of Resend Request | |
| Text Formats | 24 |
| Batch-Item Query Text | 25 |
| Item Request Text | 26 |
| Item-Selection Instruction Text | |
| Item-selection instruction text (no previous value) | |
| Item-selection instruction text (previous value exists) | |
| Measurement Data Text | |
| Mark (Result Flag) Specifications | 37 |

Purpose

This document describes the communications between the **ADVIA**[®] **1650** Chemistry System and a host computer (host).

The types of online communications between the ADVIA 1650 Chemistry System and a host are:

- The ADVIA 1650 Chemistry System requests a specific workorder or a batch of workorders from the host computer.
- The host computer sends the workorders (item selections) in reply to these requests.
- The ADVIA 1650 Chemistry System sends the sample results for patient and control samples in real time, or in batch when initiated by the operator.

The **ADVIA 1650** Chemistry System can also communicate with external samplers such as the rack handler and laboratory automation systems (**ADVIA**[®] **LabCeII**[™] System). For information about these types of communications, refer to the publication entitled <u>ADVIA 1650</u> Chemistry System and Laboratory Automation System Connection Specifications.

Interconnecting the ADVIA 1650 Chemistry System and the Host Computer

As shown below, the **ADVIA 1650** Chemistry System is connected to the host computer via an RS-232C cable supplied by the laboratory. A 25-pin (RS-232C) to 9-pin adapter is required.




Figure 1 Simplified Interconnection Diagram

Communication Specifications and Parameters

The following specifications and parameters are applicable to the interface between the ADVIA 1650 Chemistry System and the host computer. The communication parameters are entered in the On-Line Set window. The underlined parameter values are the default selections.

RS-232C Transmission specifications:

Asynchronous Synchronization:

4800 bps, 9600 bps, 14400 bps, 19200 bps Communication speed:

Coding:

Parity: None, Even, Odd

Start bit: <u>1 bit</u> Stop bit: 1 bit, <u>2 bits</u> Data length: 7 bits, 8 bits

Number of retries:

0, 1, 2, 3 Number of error text skips:

> The number of times that online processing of text alone can be skipped without the system going offline, in the event that a normal response cannot be obtained

because of an abnormality in the text format.

D-SUB 25-pin female ADVIA 1650 system connector:

(Use a male connector on the cable side.)

See Table 1 on page 3. Connector signals:

The basic procedure is alternately to transmit a message and to receive an Transmission procedure:

acknowledgment (ACK or NAK).

Online parameters

Automatic item selection (STT): No, Yes Automatic real item selection (Rack Han-No, Yes

dler or LAS):

Automatic transfer of usual sample

results:

No, Yes

Automatic transfer of interruption sample

results:

No, Yes

Automatic transfer of pilot sample results:

No, Yes

Automatic transfer of STAT sample

results:

No. Yes

Online timeout setting:

Frame interval:

The following timeout values can be set.

Automatic item selection: (Timeout time from the start to end of

automatic item selection.)

60.0 seconds

Real-time item selecting:

3.0 seconds

(Timeout time from the start to end of real-time item selection.)

3.0 seconds

(Timeout time between frames)

Determines which items (tests) are sent in the real-time output, batch output, Online item setting:

and qualitative output. You can use the ADVIA 1650 processing item numbers (1 - 326) from the **Process Sequence** window, or you can define your own item codes (host numbers) on the Item Setting dialog box of the On-Line Set window to identify the items (tests). The number that is entered into the code

field on the Item Setting Window, must be right justified.

Connector Pin Assignments

Table 1 lists the assignments for the 9-pin female connector attached to the COM1 port of the ADVIA 1650 workstation.

Table 1 9-Pin Connector Assignments

| No. | I/O | Signal Name | Connection Information | |
|-----|-----|-------------|-------------------------|--|
| 1 | IN | DCD | not connected | |
| 2 | IN | SIN | ADVIA data IN | |
| 3 | OUT | SOUT | ADVIA data OUT | |
| 4 | OUT | DTR | Connect to pin 6 (DSR). | |
| 5 | | GND | signal ground | |
| 6 | IN | DSR | Connect to pin 4 (DTR). | |
| 7 | OUT | RTS | Connect to pin 8 (CTS). | |
| 8 | IN | CTS | Connect to pin 7 (RTS). | |
| 9 | IN | RI | not connected | |

Connect pins 4 (DTR) and 6 (DSR) together, and pins 7 (RTS) and 8 (CTS) together within the 9-pin connector since the ADVIA 1650 Chemistry System itself does not control these signals.

Use a separate shielded cable. To prevent external noise from entering the system, connect the shield to the connector casing at each end of the cable (if not already done).

Host connector: Omron XM2D-2501

XM2Z-0001

9-pin adapter: Omron XM2A-2501 (soldered)

XM2S-2511 (screw-on case)

When connecting this cable to the COM1 port connector at the rear of the **ADVIA 1650** workstation (Figure 2), be sure to route it away from AC cables and any equipment that generates noise.

Do not try to connect the 25-pin RS-232C connector on the cable directly to the 9-pin COM1 connector.




Figure 2 COM1 Port Connector

NOTE

COM1 port locations vary with different models of PC:

PC Model SP6V COM1 is located on the top.
PC Model SSE350 COM1 is located on the bottom.
PC Model Seattle2 COM1 is located on the bottom.

Message Frame Format

Each message frame exchanged between the **ADVIA 1650** Chemistry System and the host computer contains the following elements:

| STX F# [0-7] Message Text ETX or ETB CSH CSL CR | LF | |
|---|----|--|
|---|----|--|

Max frame length (selectable) 256 bytes, <u>512 bytes</u>

Code used ASCII, shift JIS STX (02H) Frame start code

F# frame number "1" (31H) to "7" (37H), "0" (30H) cyclically

After ENQ transmission, the frame number is incremented beginning with 1, and when

the frame number reaches 7, it returns to 0

Message Text to be transmitted

ETX (03H) Final-frame code

ETB (17H) Intermediate-frame code

CSH Checksum upper
CSL Checksum lower
CR (0DH) + LF (0AH) Frame end code

Communication Control Codes

The following control codes provide status information:

ENQ (05H) Enquire Connection

ACK (06H) Affirmative acknowledge (acknowledgment)
NAK (15H) Negative acknowledge (acknowledgment)

EOT (04H) End of transmission

DC1 (11H) Sample request skip response (equipment control character)

The communication control codes are sent independently. (They are not enclosed by STX and ETX codes.)

Checksum Test No, Yes

The Checksum test, which is enabled or disabled using the **On-Line Set** window, checks the quality of the message. The sum of the characters, after STX, up to and including ETX or ETB, using modulo 256, is read hexadecimally and expressed as a 2-digit ASCII character. In the case of checksum **No**, a space code (20H) is transmitted, and no check is performed on the system side.

Example: <STX>1ABCDE<ETX><CSH><CSL><CR><LF>

- "1" (49D)
- "A" (65D)
- "B" (66D)
- "C" (67D)
- "D" (68D)
 "E" (69D)
- "E" (69D) · ETX (03D)

387D \rightarrow 183H (hexadecimal) \rightarrow modulo 256 \rightarrow 83H

The checksum value 83H is sent as follows: CSH = "8" (38H) and CSL = "3" (33H).

If both the sending and receiving device obtain the same checksum values, the message is considered valid.

4 Message Frame Format

Communication Modes

The available modes of communication are determined by the hardware configuration and by the specific options selected by the user.

Hardware Configuration

For communication purposes, the hardware configuration is determined by the sampling and barcode devices installed on your **ADVIA 1650** Chemistry System. This information is contained on the **System Specifications Set** window and is summarized as follows:

| Specification | Barcode Option | Sample Delivery Option |
|------------------------------|----------------------------|--|
| Sample Tray (no barcode) | Yes <u>No</u> | Sample Tray / Rack conveyer / External Transport |
| Sample Tray (with barcode) | <u>Yes</u> No | Sample Tray / Rack conveyer / External Transport |
| Rack conveyer (no barcode) | Yes <u>No</u> 1 | Sample Tray / Rack conveyer / External Transport |
| Rack conveyer (with barcode) | <u>Yes</u> No ¹ | Sample Tray / Rack conveyer / External Transport |
| Laboratory automation system | <u>Yes</u> No ² | Sample Tray / Rack conveyer / External Transport |

Table 2 System Specifications Set Hardware Configuration

System Customization

The operator can customize the following system features:

Feature Options Applicable Window Sampler for patient samples on-board sampler external sampler System Specifications Set Sample barcode reader System Specifications Set Bar-code cup posi. Automatic Item Select Yes No On-Line Set Results reporting Auto transfer, User Item Code, Data Clean On-Line Set

Table 3 System Customization

· Sampler Selection for Patient Samples

Assuming that an external sampler (Rack Handler, LAS) is configured on the **System Specifications Set** window, the **Start Conditions Set** window determines which sampler is used for patient samples in the next run. For example, if the external sampler is unavailable, the operator can use the onboard sampler (STT) to run the patient samples.

Barcode Selection

Assuming the sample barcode reader is configured on the **System Specifications Set** window, the **Start Conditions Set** window determines if the barcode mode or the position mode identifies samples on the onboard sampler.

This option is valid only for the onboard sampler. The barcode mode must be used for the external sampler.

¹ This setting is not supported. The rack handler used by the **ADVIA 1650** Chemistry System is managed as a laboratory automation device (LAS).

² Currently, all external sampling devices, including the rack handler, must have a barcode reader.

Automatic Item Select - only applicable to patient samples aspirated from the onboard sampler (STT)
 Choose Yes for the Automatic Item Select option on the On-Line Set window to have workorders for the current tray (maximum 84) downloaded automatically from the host computer at the start of the run.

Choose No when patient samples are going to be aspirated from an external sampler such as the rack handler, or when you are going to create workorders at the **ADVIA 1650** Chemistry System workstation, or when you are going to download workorders created at the host computer in batches.

Data Output

Use the **On-Line Set** window (Figure 3) to determine if sample results are sent to the host computer in real time:

Automatic transfer of general sample results (STT)

Automatic transfer of interruption sample result (Rack Handler, LAS)

Automatic transfer of pilot sample results (Controls)

Automatic transfer of urgent (STAT) sample results (STT)




Figure 3 On-Line Set Window

6 Communication Modes

Available Operating Modes

Table 4 Operating Modes

| Mode | System Configuration | Case | Sample ID | Sample feed | Order regis- tration | Main purpose |
|------|--|------|-----------------------------|------------------------------|----------------------------|---|
| 1 | Onboard sampler only (no barcode reader) | <1> | TT-CUP | Sample Tray | Auto. | Operation in which a TT-CUP number work sheet is drawn up by the host |
| 2 | | <2> | TT-CUP | Sample Tray | Batch | When automatic item registration ¹ is not possible |
| 3 | Onboard sampler only (with barcode reader) | <3> | Barcode | Sample Tray | Auto. | Sample numbers are used as barcodes |
| 4 | | <4> | Barcode | Sample Tray | Batch | When automatic item registration is not possible |
| 5 | | <1> | TT-CUP | Sample Tray | Auto. | Processing of a sample that does not permit barcode operation |
| 6 | | <2> | TT-CUP | Sample Tray | Batch | Processing of a sample that does not permit barcode operation |
| 7 | Rack handler | <5> | Barcode | Rack option | Real | Sample numbers are used as barcodes |
| 8 | (with barcode) | <6> | Barcode | Rack option | Batch | When real item registration ¹ is not possible |
| 9 | These modes are not sup- | <7> | Rack number | Rack option | Real | When barcode operation is not possible |
| 10 | ported at this time. Instead, the rack handler | <8> | Rack number | Rack option | Batch | When barcode operation is not possible |
| 11 | is managed as an LAS device as described in | <9> | SEQ number ² | Rack option | Batch | Sample numbers automatically occurred are used |
| 12 | modes 17 to 24. | <9 > | SEQ number ² | Rack option | Real | Sample numbers automatically occurred are used |
| 13 | | <3> | Barcode | Sample Tray | Auto. | Interruption sample processing |
| 14 | | <4> | Barcode | Sample Tray | Batch | Interruption sample processing |
| 15 | | <1> | TT-CUP | Sample Tray | Auto. | Processing of a sample that does not permit barcode operation |
| 16 | | <2> | TT-CUP | Sample Tray | Batch | Processing of a sample that does not permit barcode operation |
| 17 | Laboratory Automation System (LAS) device | <5> | Barcode | Laboratory automation system | Real | Sample numbers are used as barcodes |
| 18 | available including the rack handler | <6> | Barcode | Laboratory automation system | Batch | When real item registration ¹ is not possible |
| 19 | | <7> | Rack number ³ | Laboratory automation system | Real | When barcode operation is not possible |
| 20 | | <8> | Rack number ³ | Laboratory automation system | Batch | When barcode operation is not possible |
| 21 | | <3> | Barcode | Sample Tray | Auto. | Interruption sample processing |
| 22 | | <4> | Barcode | Sample Tray | Batch | Interruption sample processing |
| 23 | | <1> | TT-CUP | Sample Tray | Auto. | Processing of a sample that does not permit barcode operation |
| 24 | | <2> | TT-CUP | Sample Tray | Batch | Processing of a sample that does not permit barcode operation |

¹ Refer to page 12 for explanations of the different communication protocols.

In cases <1> and <7>, it is necessary to prepare a loadlist using the host, and also to set the data concerning the position number in the item selection instruction text.

In cases <2> and <8>, it is possible to draw up a loadlist using the ADVIA 1650 Chemistry System as well as the host.

In case (9), it is necessary to draw up a worksheet using the host, and also to set the data concerning the sample number in the item selection instruction text.

² Sequential No. analysis

³ This setting is not presently supported.

Sample Identification

The **ADVIA 1650** Chemistry System has the following means of identifying samples internally within its own database, during sampling, and during the transfer of data between the system and the host computer or an external sampler (rack handler or other LAS device).

Registration Number

The registration number is a sequence number between 1 and 10000. This number is employed in various **ADVIA 1650** workstation windows, but it cannot be used to identify samples in the communications between the **ADVIA 1650** Chemistry System and the host computer or an external sampler.

Position Number

When not using barcodes, the sample position number is used to identify patient samples aspirated from the onboard sampler (STT). You can also use the position number as a data inquiry key in various **ADVIA 1650** workstation windows.

The STT position number format is: 01-01 TT-CC (TT= tray number, 1 to 97; CC = sample position, 1to 84)

Position numbers are used in online communications messages. The numbers are left-justified, not zero-suppressed, and terminated with an ASCII Space character (20H).

When position numbers are used, the entry of a sample number is still required.

Position numbers cannot be used to identify samples aspirated from an external sampler such the rack handler.

Sample Number

The sample number consists of up to 13 alphanumeric characters. When barcode labels are used, the sample number becomes the barcode ID that identifies patient samples aspirated from the onboard sampler (STT) or an external sampler. You can also use the sample number as a data inquiry key in various **ADVIA 1650** workstation windows.

A sample number is required for all workorders. This includes the position mode as described in the previous paragraph. The **ADVIA 1650** Chemistry System will not accept a workorder that is missing a sample number.

The sample number formats are:

- General samples and interruption samples
 - You can use IDs consisting of up to 13 arbitrary characters.
- For STAT samples

A sample number from E01 to E84 (3 characters, fixed length) is automatically assigned when the operator uses the STAT button on the Operation Panel to run a STAT sample.

Control samples

A sample number from PA01 to PJ20 (4 characters, fixed length) is automatically assigned when a control sample is run.

Pxnn (P = pilot [control], x = control letter assigned on QC Sample Definition window, nn = sequence number)

Sample numbers are used in online communications messages. When transmitted, these sample numbers are left-justified, not zero-suppressed, and terminated with an ASCII Space character (20H). For example: "001" and "1" are recognized as different sample numbers.

Sample Information

The sample information consists of comment 1, comment 2, sex, age, and so on.

Comments 1 and 2 each consist of a maximum of 16 arbitrary characters. They can be used as patient numbers or patient names. They can also be used as search keys for making data inquiries. There is no check for duplicate entries.

In online communications messages, the sample data is treated as character data.

8 Communication Modes

Classification of ADVIA 1650 Samples

Table 5 Sample Types

| Sample | | Item | Data Output | Sample setting ¹ | |
|----------------------------------|------------------------------------|--------------------------|-------------|-----------------------------|---------------------|
| | Sample ID | Registration by the Host | to the Host | Sample Tray | External Sampler |
| General sample ² | 13 digits | yes ³ | yes | yes ⁴ | yes ⁵ |
| Interruption sample ⁶ | 13 digits | yes ³ | yes | yes | no |
| STAT sample ⁷ | E01 – E84 | no | yes | yes | no |
| Control sample ⁸ | PA01 – PJ20 | no | yes | yes | no |
| Calibration ⁹ | C0101 – C0110 M980101 – M980110 | no | no | yes | no |

All samples except general samples must be run on the Sample Tray (onboard sampler).
General samples are run on the Sample Tray or an external sampler such as the rack handler or a laboratory automation system (LAS), if available. When the run is started, the operator uses the **Start Conditions** window to select which sampler is used for general samples. If for any reason multiple samplers are selected, the analysis priority is to use the Sample Tray first and then use the rack handler or LAS.

Single-point calibration: Cxxyy xx: Cup No. (01-61) yy: Number of times (01-05)

Multi standard: Mxxyyzz xx: TT No. (98,99) yy: Cup No. (01-84) zz: Number of times (01-05)

² The **ADVIA 1650** Chemistry System can manage a maximum of 10,000 general samples.

³ The host computer can create workorders (item registration) for general samples and interruption samples.

⁴ The position number for a sample on the Sample Tray (STT) consists of a tray number and a cup number. The tray number (TT) must be from 01 to 97, while the cup number (CUP) must be from 01 to 84.

⁵ Currently, rack position numbers are not used for external samplers. Barcode labels are required for the LAS devices including the rack handler. However for reference purposes, the rack position number consists of a rack number and a sample number. The rack number must be from 0001 to 9999, while the sample number must be from 01 to 10.

⁶ Interruption samples are treated as general samples; however, these samples are set at the CUP positions preset on the Sample Tray, and are analyzed first.

⁷ STAT samples are also set at predetermined positions on the Sample Tray, and are analyzed preferentially. However, STAT item selection cannot be downloaded from the host computer. STAT test selectivity is determined using the **STAT Order Setup** window. The sample IDs E01 to E84 are reserved for STAT samples.

⁸ Quality control samples are aspirated from positions on the Sample Tray (CTT). The listed sample IDs are reserved for controls only. Controls cannot be requested from the host computer. Instead, controls are usually requested from the **Start Conditions** window.

⁹ Calibration samples are not handled online.

Item Order Operation

When the host computer is online, all of the sample orders, regardless of analysis results of the sample, are inquired by the host computer. Then, workorders from the host computer and the **ADVIA 1650** Chemistry System are combined, and the **ADVIA 1650** Chemistry System selects the analysis orders.

The following Table explains how the system manages workorders from the host computer.

Table 6 Managing Workorders from the Host Computer

| Host Order | System Order | System Response |
|--------------------------|--------------|---|
| New workorder | No | The system registers the sample and runs the tests requested by the host computer as a new workorder. |
| New workorder | Yes | The sample in the system is deleted, and the tests requested by the host computer are run as a new workorder. |
| Addition, Re-analysis | No | The sample is registered in the system, and the tests requested by the host computer are run as a new workorder. |
| Addition, Re-analysis | Yes | If tests requested by the host computer have already been run, they are repeated using the re-analysis conditions. |
| | | If tests requested by the host computer have not been run, they are run as originals under first-analysis conditions. |
| No workorder | No | This sample is not analyzed. |
| No workorder | Yes | The tests in the system-created workorder are run. |
| Sample deletion | No | This sample is not analyzed. |
| Sample deletion | Yes | System order information of this sample is deleted. No analysis. |

10 Communication Modes

Host Item Code Table

The ADVIA 1650 Chemistry System can manage a maximum of 326 items:

- 300 test items
- 20 computed result items (ratios)
- 3 serum indices (lipemia, hemolysis and icterus)
- 3 electrolyte tests (Na⁺, K⁺, and Cl⁺)

There are two ways to identify these tests for online communications. You can use the processing item numbers employed by the **ADVIA 1650** Chemistry System, or you can define a separate set of item numbers or host item codes that agree with your laboratory's numbering scheme.

To enter the user-defined item codes, first select **Do** for the Used Item Code option on the **On-Line Set** window, and then enter the required numbers in the Code boxes of the Item Setting dialog box. Thereafter, you can use the host item codes when ordering tests or managing results at the host computer. However, you must use the system processing numbers when working at the **ADVIA 1650** Chemistry System.

The Serum info.item area in the **System Specification Set** window determines if the serum indices are automatically run for each sample or if they must be requested each time using the **Order Entry** window.

When any electrolyte test is ordered, the other two are run automatically.

In the indication of the result of each item, the number of digits to the right the decimal point can be adjusted for each item. Select the number of digits.

Table 7 Host Item Codes

| Processing Item No. | Host Item Codes | Item Name | Number of Digits (Indication) | To right of decimal point | Classification |
|------------------------|--------------------|-------------|-------------------------------|---------------------------|------------------|
| xx1 | xx1 ¹ | (AST, etc.) | 8 digits | test specific | photometric test |
| xx2 | xx2 | " | 8 digits | test specific | photometric test |
| xxx | xx3 | " | 8 digits | test specific | photometric test |
| xxx | xxx | " | 8 digits | test specific | photometric test |
| xxx | xxx | " | 8 digits | test specific | photometric test |
| xxx | xxx | " | 8 digits | test specific | photometric test |
| xxx | XXX | A/G, etc. | 8 digits | test specific | ratio |
| xxx | xxx | II | 8 digits | test specific | ratio |
| xxx | XXX | " | 8 digits | test specific | ratio |
| xxx | XXX | " | 8 digits | test specific | ratio |
| xxx | xxx | II | 8 digits | test specific | ratio |
| xxx | XXX | Lipemia | Qualitative display possible | test specific | Blood serum data |
| XXX | xxx | Hemolysis | Qualitative display possible | test specific | Blood serum data |
| XXX | xxx | Icterus | Qualitative display possible | test specific | Blood serum data |
| XXX | xxx | Na | 8 digits | test specific | Electrolyte test |
| XXX | xxx | K | 8 digits | test specific | Electrolyte test |
| XXX | xxx | CI | 8 digits | test specific | Electrolyte test |

¹ Normal item code is zero-suppressed, right-justified.

Communication Protocols

The following communication protocols are available:

Batch-item registration

Communication starts after the operator uses the Host Request dialog box of the **Order Entry** window to identify a group (batch) of workorders to be downloaded from the host computer.

Automatic-item registration (Host Query Mode)

When this feature is selected on the **On-Line Select** window, the host computer downloads available workorders for STT patient samples upon request at the start of the run.

· Real-item registration

This communications protocol supports external samplers such as the rack handler or a laboratory automation system. The host computer downloads an available workorder for each sample identification number supplied by the external sampler.

Batch-data output

Communication starts after the operator uses the Host Transfer dialog box of the **Review / Edit** window to identify a group (batch) of sample results to be sent to the host computer.

Real-data output

Whenever each sample's assay is completed, the results are sent to the host computer.

The following information is applicable to each of the communication protocols.

Time outs

The system-side timer value is the maximum time that the **ADVIA 1650** Chemistry Systems waits for a correct response from the host. Similarly, the host-side timer is the maximum time allowed for the **ADVIA 1650** Chemistry System to respond.



IMPORTANT

The **ADVIA 1650** Chemistry System can analyze individual items at three-second intervals. If the sum of the system-side and host-side values exceeds the analysis time interval, the system sampling time is increased resulting in decreased throughput. To avoid a significant decrease in throughput, the host-side timer value should be less than that of the analyzer, and the host computer should respond as quickly as possible.

If there are 'n' request items for one sample, the sample analysis interval is 3.0 seconds x 'n'. Consequently, the greater the value of n, the smaller the possibility of decreased throughput during the communication waiting period.

The three-second cycle time of the **ADVIA 1650** imposes some restrictions on the response time of the Host computer program. If one wants to continuously process samples with a single photometric test requested on the work order (single item registration) then the Host computer response time must be less than 3.0 seconds. A response time greater than 3.0 seconds will reduce system throughput. The **ADVIA 1650** must wait for a Host response before sampling can continue.

In addition the **ADVIA 1650** gives priority to the communication of work orders (item registration) over result transfer (data communication). If the item registration is not completed within 3.0 seconds then data will accumulate in the **ADVIA 1650** PC. This accumulation of data can damage the ADVIA PC.

The need for a 3.0 second response is reduced if two or more photometric tests are requested per work order. The Host reponse time should be less than 3.0 seconds times the number of tests requested (3.0 sec x 'n' tests). ISE analysis by the **ADVIA 1650** requires 24 seconds.

Table 8 Example of Response Time to Meet 3.0 Second Limitation

| On condition of t | he transmissi | on rate 9600bps | | | |
|--------------------------|-----------------|-----------------------|-----------------------------|---|---|
| | | Communication Byte | Communciation Time (sec) | ADVIA Typical Response Time (sec) | Host Estimate Response Time (sec) |
| Item Registration | | | | | |
| $ADVIA \rightarrow HOST$ | ENQ | 1 | 0.00104 | | |
| $ADVIA \leftarrow HOST$ | ACK | 1 | 0.00104 | | 0.25 |
| $ADVIA \rightarrow HOST$ | Item request | 30 | 0.03125 | 0.13 | |
| ADVIA ←HOST | ACK | 1 | 0.00104 | | 0.25 |
| $ADVIA \rightarrow HOST$ | EOT | 1 | 0.00104 | 0.13 | |
| $ADVIA \leftarrow HOST$ | ENQ | 1 | 0.00104 | | 0.25 |
| $ADVIA \rightarrow HOST$ | ACK | 1 | 0.00104 | 0.13 | |
| $ADVIA \leftarrow HOST$ | Item regist. | 93 | 0.09687 | | 0.25 |
| $ADVIA \rightarrow HOST$ | ACK | 1 | 0.00104 | 0.13 | |
| ADVIA ← HOST | EOT | 1 | 0.00104 | | 0.25 |
| Data Online | | | | | |
| $ADVIA \rightarrow HOST$ | ENQ | 1 | 0.00104 | | |
| ADVIA ← HOST | ACK | 1 | 0.00104 | | 0.25 |
| $ADVIA \rightarrow HOST$ | Data online | 112 | 0.11667 | 0.13 | |
| ADVIA ← HOST | ACK | 1 | 0.00104 | | 0.25 |
| $ADVIA \to HOST$ | ЕОТ | 1 | 0.00104 | 0.13 | |
| | | 247 | 0.25729 | 0.78 | 1.75 |
| 7 | Total Online Ti | me | | 2.79 | |

Required host response to communication initialization

The **ADVIA 1650** Chemistry System initiates communications regardless of the receiving condition of the host computer. The first line request (ENQ) of the communication procedure is communicated to identify whether the host computer can receive or not. Regardless of the host process conditions, the host must receive the line request (ENQ) and reply for confirmation.

If the host computer is temporarily unable to respond, the **ADVIA 1650** Chemistry System will retransmit a line request (ENQ) repeatedly up to the maximum number of times¹, if a negative response (NAK) is returned to the **ADVIA 1650** Chemistry System line request within an analyzer-side time-out. This retransmission can expand the host reception starting time. Note that these processes reduce system performance.

Host means of identifying samples

The position-mode method (TT-CUP No.) is available for patient samples aspirated from the on-board sampler (STT). For example, a user who operates the **ADVIA 1650** Chemistry System in the barcode mode can request a suitable ID that is different from the barcode ID during analysis, and process interruption samples using the TT-CUP numbers. If, however, the barcode sample is still being analyzed, the **ADVIA 1650** Chemistry System cannot be offline. In this case, the automatic item registration by TT-CUP No. message is transmitted to the host. The ID that the host does not handle is transmitted, and the host must respond as follows:

To item inquiries: No order response Confirmation response

If the response is not transmitted, a time-out error occurs on the analyzer side.

¹ Refer to the error processing information on page 22.

Batch Item Registration Protocol

Communication is started from the **Order Entry** window.

The request conditions (sample identification, range and related data) entered by the operator are inserted into the batch item inquiry text which is then transferred to the host computer as a workorder inquiry.

The host sends the item-selection instruction text, under specified conditions, for the samples with items to be analyzed by the **ADVIA 1650** Chemistry System.

If these communications are interrupted by the operator, the **ADVIA 1650** Chemistry System sends an EOT signal to the host, and the transmission is ended.

The event numbers in the following diagram, for example <1>, are explained in Table 9.

| System | ENQ | \rightarrow | | Host |
|--------|--------------------|---------------|---|------|
| <1> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | Batch Item inquiry | \rightarrow | | Host |
| <3> | | | | |
| System | | ← | Acknowledgment | Host |
| <4> | | | | |
| System | EOT | \rightarrow | | Host |
| <5> | | | | |
| System | | ← | ENQ | Host |
| <6> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| <7> | | | | |
| System | | ← | Item-selection instruction ¹ | Host |
| <8> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| | | (Repetition) | | |
| System | | ← | Item-selection instruction | Host |
| <8> | | | | |
| System | Acknowledgment | → | | Host |
| <7> | | | | |
| | | ← | EOT | Host |

Figure 4 Batch Item Registration Dialog

¹ Only samples that have item-selection instructions are registered. If the line was finished on this first timing, there is no registration.

Table 9 Batch Item Registration Events

| | Timer | Timer Monit | oring Interval | Timer Interval | Action |
|-----|--|-------------------------|----------------------------------|---------------------------------------|---|
| | imer | Start | End | - Timer interval | Action |
| <1> | ENQ | ENQ sent | Receipt of acknowledgment | 5 seconds | EOT is sent, and processing stops. |
| <2> | System-side processing period | Acknowledgment received | Inquiry sent | (2 seconds) | |
| <3> | Inquiry acknowledg- ment | Inquiry sent | Acknowledgment received | 3 seconds | EOT is sent, and processing stops. |
| <4> | System-side processing period | Acknowledgment received | EOT sent | (2 seconds) | |
| <5> | Response after EOT | EOT sent | ENQ received | 5 seconds + a ¹ changeable | The system sends EOT, and item selection processing stops. |
| <6> | System-side processing period | ENQ received | Acknowledgment sent | (2 seconds) | |
| <7> | Response after an acknowledg- ment | Acknowledgment sent | Next-frame mes- sage received | 5 seconds + a changeable | The system sends EOT, and item selection processing is interrupted. The received part is valid. |
| <8> | System-side processing period | Item selection received | Acknowledgment sent | 2 seconds | |
| <9> | Entire procedure processing time | System ENQ sent | Host EOT sent | Not observed | |

¹ a is the time-out value between frames.

Automatic Item Registration Protocol

If the Automatic item select option is selected on the **On-line Set** window, the system requests workorders for the current tray (samples loaded onto the STT). When the Sample Tray-analysis mode starts, communications take place and analysis starts automatically after the workorders are downloaded and registered.

The event numbers in the following diagram, for example <1>, are explained in Table 10.

| System | ENQ | \rightarrow | | Host |
|--------|----------------|---------------|---|------|
| <1> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | Item inquiry | \rightarrow | | Host |
| <3> | | | | |
| System | | ← | Acknowledgment | Host |
| <4> | | | | |
| System | EOT | \rightarrow | | Host |
| <5> | | | | |
| System | | ← | ENQ | Host |
| <6> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| <7> | | | | |
| System | | ← | Item-selection instruction ¹ | Host |
| <8> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| 1 | | (Repetition) | | 1 |
| System | | ← | Item-selection instruction | Host |
| <8> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| <7> | | | | |
| | | ← | EOT | Host |

Figure 5 Automatic Item Registration Dialog

Only samples that have item-selection instructions are registered. If the line was finished on this first timing, there is no registration.

Table 10 Automatic Item Registration Events

| | Timer Monitoring Inter | | oring Interval | - Timer Interval | Action |
|-----|--|-------------------------|----------------------------------|--|------------------------------------|
| | IIIIei | Start | End | - Tilliel lillerval | Action |
| <1> | ENQ | ENQ sent | Receipt of acknowledgment | 5 seconds | EOT is sent, and processing stops. |
| <2> | System-side processing period | Acknowledgment received | Inquiry sent | (2 seconds) | |
| <3> | Inquiry acknowledg- ment | Inquiry sent | Acknowledgment received | 3 seconds | EOT is sent, and processing stops. |
| <4> | System-side processing period | Acknowledgment received | EOT sent | (2 seconds) | |
| <5> | Response after EOT | EOT sent | ENQ received | 5 seconds + a ¹ changeable | EOT is sent, and processing stops. |
| <6> | System-side processing period | ENQ received | Acknowledgment sent | (2 seconds) | |
| <7> | Response after an acknowledg- ment | Acknowledgment sent | Next-frame mes- sage received | 5 seconds + a changeable | EOT is sent, and processing stops. |
| <8> | System-side processing period | Item selection received | Acknowledgment sent | 2 seconds | |
| <9> | Entire procedure processing time | System ENQ sent | Host EOT sent | 42 seconds + b ² changeable | EOT is sent, and processing stops. |

¹ a is the time-out value between frames.

² Automatic item-selection time-out value

Real Item Registration Protocol

During analysis in which patient samples are aspirated from an external sampler such as the rack handler or LAS device, the external sampler reads the barcode ID for each sample and sends it to the **ADVIA 1650** Chemistry System. If the required workorder is not available, the system queries the host computer.

The event numbers in the following diagram, for example <1>, are explained in Table 11.

| System | ENQ | → | | Host |
|--------|----------------|---------------|---|----------|
| <1> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | Item inquiry | → | | Host |
| <3> | | | | |
| System | | ← | Acknowledgment | Host |
| <4> | | | | |
| System | EOT | → | | Host |
| <5> | | | | |
| System | | ← | ENQ | Host |
| <6> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| <7> | | | | |
| System | | ← | Item-selection instruction ¹ | Host |
| <8> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| - | | (Repetition) | | . |
| System | | ← | Item-selection instruction | Host |
| <8> | | | | |
| System | Acknowledgment | \rightarrow | | Host |
| <7> | | | | |
| | | ← | EOT | Host |

Figure 6 Real Item Registration Dialog

Only samples that have item-selection instructions are registered. If the line was finished on this first timing, there is no registration.

Table 11 Real Item Registration Events

| | Timer Monitoring Inter | | oring Interval | - Timer Interval | Action |
|-----|--|-------------------------|----------------------------------|--|------------------------------------|
| | IIIIei | Start | End | - Tilliel lillerval | Action |
| <1> | ENQ | ENQ sent | Receipt of acknowledgment | 5 seconds | EOT is sent, and processing stops. |
| <2> | System-side processing period | Acknowledgment received | Inquiry sent | (2 seconds) | |
| <3> | Inquiry acknowledg- ment | Inquiry sent | Acknowledgment received | 3 seconds | EOT is sent, and processing stops. |
| <4> | System-side processing period | Acknowledgment received | EOT sent | (2 seconds) | |
| <5> | Response after EOT | EOT sent | ENQ received | 5 seconds + a ¹ changeable | EOT is sent, and processing stops. |
| <6> | System-side processing period | ENQ received | Acknowledgment sent | (2 seconds) | |
| <7> | Response after an acknowledg- ment | Acknowledgment sent | Next-frame mes- sage received | 5 seconds + a changeable | EOT is sent, and processing stops. |
| <8> | System-side processing period | Item selection received | Acknowledgment sent | 2 seconds | |
| <9> | Entire procedure processing time | System ENQ sent | Host EOT sent | 42 seconds + b ² changeable | EOT is sent, and processing stops. |

¹ a is the time-out value between frames.

² Automatic item-selection time-out value

Batch Data Output Protocol

Communication is started from the **Review/Edit** window, where the operator identifies the group of sample results to be sent to the host computer.

You can use the Data Clean Set feature of the **On-Line Set** window to validate the sample results before they are transmitted. Results that fail are not sent. Use the Test Results List dialog box of the **Review / Edit** window to determine which samples failed the data clean check.

If the operator interrupts this communication, the **ADVIA 1650** Chemistry System sends an EOT signal to the host and the transmission is ended.

The event numbers in the following diagram, for example <1>, are explained in Table 12.

| System | ENQ | \rightarrow | | Host |
|--------|----------------|---------------|----------------|----------|
| <1> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | Sample results | \rightarrow | | Host |
| <3> | | | | |
| System | | ← | Acknowledgment | Host |
| | | (Repetitio | n) | <u> </u> |
| System | Sample results | \rightarrow | | Host |
| <3> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | EOT | \rightarrow | | Host |

Figure 7 Batch Data Output Dialog

Table 12 Batch Data Output Events

| | Timer | Timer Monitoring Interval | | Timer Interval | Action |
|-----|-------------------------------|---------------------------|---------------------------|----------------|------------------------------------|
| | Timer | Start | | | Action |
| <1> | ENQ | ENQ sent | Receipt of acknowledgment | 5 seconds | EOT is sent, and processing stops. |
| <2> | System-side processing period | Acknowledgment received | Sample results sent | (3 seconds) | |
| <3> | Data acknowl- edgment | Sample results sent | Acknowledgment received | 3 seconds | EOT is sent, and processing stops. |
| <4> | All processing period | | | not observed | |

Real Data Output Protocol

The results for each sample are automatically sent to the host computer when they are available.

Use the Automatic transfer area on the **On-Line Set** window to designate which kinds of sample results are sent. Use the Data Clean Set feature to validate the sample results before they are transmitted. Results that fail are not sent. Use the Test Results List dialog box of the **Review / Edit** window to determine which samples failed the data clean check.

The event numbers in the following diagram for example <1>, are explained in Table 13.

| System | ENQ | \rightarrow | | Host |
|--------|----------------|---------------|----------------|------|
| <1> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | Sample results | \rightarrow | | Host |
| <3> | | | | |
| System | | ← | Acknowledgment | Host |
| <2> | | | | |
| System | EOT | \rightarrow | | Host |

Figure 8 Real Data Output Dialog

Table 13 Real Data Output Events

| | Timer | Timer Monito | oring Interval | Timer Interval | Action | | |
|-----|-------------------------------|-------------------------|---------------------------|----------------|------------------------------------|------------------|--------|
| | imei | Start End | | | | Tiller lillervar | Action |
| <1> | ENQ | ENQ sent | Receipt of acknowledgment | 5 seconds | EOT is sent, and processing stops. | | |
| <2> | System-side processing period | Acknowledgment received | Sample results sent | 3 seconds | | | |
| <3> | Data acknowl- edgment | Sample results sent | Acknowledgment received | 3 seconds | EOT is sent, and processing stops. | | |
| <4> | All processing period | | | not observed | | | |

Error Processing

Abnormal Text Format

If the system receives text that is abnormally formatted, the ADVIA 1650 Chemistry System sends a NAK signal.

The abnormal text-NAK cycle may repeat up to three times.

When the system receives text that is abnormally formatted again, the error-skip counts (by the internal counter) increase, and the number of counts is compared to the set error-sample skip number. Then, either of the following processes takes place (refer to example below).

• If the number of error-skip counts is less than or equal to the number of text error-skip counts:

System Transmits request-skip response (DC1).

Host Transmits item-selection instruction for next sample, or EOT.

• If the number of error-skip counts exceeds the number of text error-skip counts:

System Transmits EOT and Item-selection processing is terminated. Valid workorders are implemented.

Host Transmission of item-selection instruction is terminated.

[Example]

Responses from the **ADVIA 1650** Chemistry System (when text that has an abnormal format is received while executing an item-selection instruction).

| | ← | Item selection-instruction | Host |
|------------------------------------|---|--|---|
| Acknowledgment (ACK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| Negative acknowledge (NAK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| Negative acknowledge (NAK) | \rightarrow | _ | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| Negative acknowledge (NAK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| (The number of error-skip counts d | oes not exc | eed the number of text error-skip counts.) | |
| Sample request skip response (DC1) | \rightarrow | | Host |
| | ← | Next sample item selection-instruction | Host |
| | • • • | | |
| | ← | Item selection-instruction | Host |
| Acknowledgment (ACK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| Negative acknowledge (NAK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| Negative acknowledge (NAK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| Negative acknowledge (NAK) | \rightarrow | | Host |
| Text format abnormal | ← | Item selection-instruction | Host |
| (The number of error-skip cour | nts exceeds | the number of text error-skip counts.) | |
| Circuit interruption (EOT) | \rightarrow | | Host |
| | Text format abnormal Negative acknowledge (NAK) Text format abnormal Negative acknowledge (NAK) Text format abnormal Negative acknowledge (NAK) Text format abnormal (The number of error-skip counts of sample request skip response (DC1) Acknowledgment (ACK) Text format abnormal Negative acknowledge (NAK) Text format abnormal (The number of error-skip cour | Acknowledgment (ACK) Text format abnormal Negative acknowledge (NAK) Text format abnormal Negative acknowledge (NAK) Text format abnormal Negative acknowledge (NAK) Text format abnormal ← (The number of error-skip counts does not exc Sample request skip response (DC1) ← Acknowledgment (ACK) Text format abnormal ← Negative acknowledge (NAK) Text format abnormal ← (The number of error-skip counts exceeds | Acknowledgment (ACK) → Item selection-instruction Negative acknowledge (NAK) → Item selection-instruction (The number of error-skip counts does not exceed the number of text error-skip counts.) Sample request skip response (DC1) → Next sample item selection-instruction ∴ Next sample item selection-instruction Acknowledgment (ACK) → Item selection-instruction Negative acknowledge (NAK) → Item selection-instruction |

Receipt Of Resend Request

If the **ADVIA 1650** Chemistry System receives a negative acknowledgment (NAK) from the host, it resends the same frame three times. At the fourth negative acknowledgment (NAK), the **ADVIA 1650** Chemistry System sends EOT and terminates the processing.

Text Formats

The following text formats specify how workorder and result information will be transmitted.

- · batch-item query text
- item-query text
- item-selection instruction text (previous time No/previous time Yes)
- measurement-data text

Typically, this information will require more than one frame:

termination code for intermediate frame = ETB (17H) + checksum + CR + LF

termination code for last frame = ETX (03H) + checksum + CR + LF

Batch-Item Query Text

(from ADVIA 1650 Order Entry window)

Applicable communication function: Batch-item registration (page 14)

Communication direction: ADVIA 1650 Chemistry System → Host Computer

Table 14 Batch-Item Query Text (fixed length, 41 bytes)

| Data Item | Bytes | Code |
|---------------------------------|-------|---|
| | | |
| Starting code | 1 | STX |
| Frame number | 1 | "1" – "7", "0" |
| Text classification | 1 | 'q' (71H) |
| Equipment identification number | 1 | Not used, space (20H) |
| Total number of blocks | 2 | "01" (no zero suppression) |
| Block number | 2 | "01" (no zero suppression) |
| ID classification | 1 | "0" (30H) Arbitrary sample ID, Barcode ID, Sequential ID "1" (31H) Sample Tray TT No. and cup position No "2" (32H) Rack No. and number indicating position in rack "3" (33H) Independent holder No |
| Initial sample ID | 13 | The response to the ID classification (previous item) is as follows: |
| | | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified, 20H at the end |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |

Table 14 Batch-Item Query Text (fixed length, 41 bytes) (Cont)

| Data Item | Bytes | Code |
|------------------|-------|---|
| Final sample ID | 13 | The response to the ID classification (previous item) is as follows: Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end Rack number: Fixed-length, 7 digits, no zero suppression, left-justified, 20H at the end Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |
| Spare | 1 | Not used (20H) |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF |

Before communication is initiated, the following information is provided by the operator using the Host request dialog box of the **Order Entry** window:

- The ID classification item indicates the method used to identify the starting sample ID and the final sample.
- The starting sample ID and the final sample ID determine the range of workorders being requested.

The **ADVIA 1650** Chemistry System does not check for the order ID and the number of received samples, relative to the item-selection instruction text of the host.

In the case of transmission without range designation, the ID positions for the initial and final samples are blank spaces.

When inspection without range designation is performed, both the blank spaces and zeros should be taken into account.

Item Request Text

Applicable communication function: Automatic item registration (page 16)

Real item registration (page 18)

Communication direction: ADVIA 1650 Chemistry System → Host Computer

| Maximum frame size | Maximum number of blocks |
|--------------------|--------------------------|
| 256 bytes | 5 blocks |
| 512 bytes | 3 blocks |

Table 15 Item Request Text (first block to last block, variable length)

| Data Item | Bytes | Code |
|--------------------------------------|--------------------------|---|
| Starting code | 1 | STX |
| Frame number | 1 | "1" – "7", "0" |
| Text classification | 1 | "Q" (51H) |
| Equipment identification number | 1 | Not used, space (20H) |
| Total number of blocks | 2 | "01" to "05" (no zero suppression) |
| Block number | 2 | "01" to "05" (no zero suppression) |
| Number of inquiry samples in a block | 2 | "01" to "99" (no zero suppression) |
| ID classification | 1 | "0" (30H) Arbitrary sample ID, Barcode ID, Sequential ID "1" (31H) Sample Tray TT No. and cup position No. "2" (32H) Rack No. and number indicating position in rack "3" (33H) Independent holder No. |
| Sample ID | 13 byte x n ¹ | The response to the ID classification (previous item) is as follows: |
| | | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified, 20H at the end |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |
| Spare | 1 | Not used (20H) |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF |

When the maximum frame size is 256 bytes, the maximum number of items is n = 18. The maximum number of samples that can be inquired with one text is 84 samples (Sample Tray) and 5 or 10 samples (rack).

When the maximum frame size is 512 bytes, the maximum number of items is n = 38. There is no item inquiry for STAT samples or quality control samples.

■ Item-Selection Instruction Text

There are two kinds of item-selection instruction text, one with previous data and the other without previous data.

The number of items for which a request can be made for one sample is 326.

Item-selection instruction text – no previous value

Applicable communication function: Batch-item registration (page 14)

Automatic item registration (page 16) Real item registration (page 18)

Communication direction: ADVIA 1650 Chemistry System ← Host Computer

| Maximum Frame Size | Maximum number of blocks |
|-----------------------|--------------------------|
| 256 bytes | 7 blocks |
| 512 bytes | 3 blocks |

Table 16 Item-selection Text - No Previous Value (first block, variable length)

| Data Item | Bytes | Code |
|---------------------------------|-------|--|
| Starting code | 1 | STX |
| Frame number | 1 | "1" – "7", "0" |
| Text classification | 1 | "O" (4FH) |
| Equipment identification number | 1 | Not used, space (20H) |
| Total number of blocks | 2 | "01" to "07" (no zero suppression) |
| Block number | 2 | "01" (no zero suppression) |
| Number of items in a block | 3 | "001" to "999" (no zero suppression) |
| Sample classification | 1 | "N" (4EH) General sample "I" (49H) Interruption sample |
| Registration data | 1 | "0" (30H) New request. (If a registered sample [workorder] already exists, it is overwritten.) "1" (31H) Item addition, re-run. (In case of an unregistered sample, same as a new request.) "2" (32H) No request. (The system does not perform registration.) "3" (33H) Sample deletion. (Deletion of a previously registered sample) |
| Sample ID | 13 | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end. Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end. Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end. Be sure to set unique sample numbers, even for an order of position No. designation. |

Table 16 Item-selection Text - No Previous Value (first block, variable length) (Cont)

| Data Item | Bytes | Code |
|--------------------------|----------------------|---|
| Position number | 7 | Filled with blank spaces when not used (20H). |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |
| Comment 1 | 16 | ASCII or shift JIS. Filled with blank spaces (20H) when not used. |
| Comment 2 | 16 | ASCII or shift JIS. Filled with blank spaces (20H) when not used. |
| Sex | 1 | "M" (4DH) male or "F" (46H) female: When not used, set to "M". |
| Age | 3 | 999 (right-justified): When not used, all three characters are 20H. Zero years old is 20H, 20H, "0". |
| Blood sampling date | 8 | YYYYMMDD Example: 19990229 Filled with blank spaces (20H) when not used. |
| Dilution coefficient | 4 | Format: 99.9 (right-justified): When not used, be sure to insert Format: ^1.0 (^0.1 - 99.9). ^ : blank space |
| Sample classification | 1 | "1" (31H) blood serum, "2" (32H) urine: When not used, be sure to input "1". |
| Container classification | 1 | "1" to "9" : When not used, "1". |
| Request item | (3+1)xn ¹ | Item number 3 bytes, Format: 999 (right-justified), 0 suppress |
| | | Sample volume and dilution determined via the Analytical Parameters (Chemistry) window. One byte "M" or "D" or "U". Normally set to "M". |
| | | "M" = Analytical conditions values |
| | | "D" or "U" = Reanalysis condition values |
| Spare | 1 | Not used (20H) |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF |

When the maximum frame size is 256 bytes the maximum number of items is n = 41. When the maximum frame size is 512 bytes the maximum number of items is n = 105.

Table 17 Item-selection Text – No Previous Value (second to last block, variable length)

| Data Item | Bytes | Code |
|---------------------------------|-------|--------------------------------------|
| Starting code | 1 | STX |
| Frame number | 1 | "1" – "7", "0" |
| Text classification | 1 | "O" (4FH) |
| Equipment identification number | 1 | Not used, space (20H) |
| Total number of blocks | 2 | "02" to "07" (no zero suppression) |
| Block number | 2 | "02" to "07" (no zero suppression) |
| Number of items in a block | 3 | "001" to "999" (no zero suppression) |
| Spare | 2 | Not used (20H) |

Table 17 Item-selection Text - No Previous Value (second to last block, variable length) (Cont)

| Data Item | Bytes | Code |
|------------------|----------------------|--|
| Sample ID | 13 | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end |
| | | Be sure to set unique sample numbers, even for an order of position No. designation. |
| Position number | 7 | Filled with blank spaces when not used (20H). |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |
| Request item | (3+1)xn ¹ | Item number 3 bytes, Format: 999 (right-justified), 0 suppress |
| | | Sample volume and dilution determined via the Analytical Parameters (C.hemistry) window. One byte "M" or "D" or "U". Normally set to "M." |
| | | "M" = Analytical conditions values |
| | | "D" or "U" = Reanalysis condition values |
| Spare | 1 | Not used (20H) |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF |

When the maximum frame size is 256 bytes, the maximum number of items is n = 54. When the maximum frame size is 512 bytes, the maximum number of items is n = 118.

Operational Notes:

- If the ID is classified as a sample ID for a batch-item inquiry or an item inquiry, set the sample ID in the itemselection instruction text, and leave the position number unused (barcode operation or creation of worksheets using the **ADVIA 1650** Chemistry System).
- If the ID is classified as a Sample Tray ID or rack-number ID for a batch-item inquiry or an item inquiry, set the position number to which the inquiry was addressed, and also set the sample number to be controlled by the host in the sample ID item (creation of worksheet using the host).
- At the start of Sample Tray analysis, the sample classified "I" as an interruption sample on sample classification is analyzed first. The ADVIA 1650 Chemistry System also has a built-in function to give priority to interruption samples for Sample Tray analysis. When the samples are analyzed on the rack option or conveyor-belt system, analysis begins sequentially as ordered regardless of that function. Therefore, interruption samples for these systems should be set on the Sample Tray.
- Registration data box

"0" New Request

A new request is unconditionally registered as a new item selection. In case of an already registered sample, its information is deleted and overwritten, thereby erasing the analysis-result data. If you cannot judge whether it is a new request or registered item on the host side, specify all of the items in the next box as "1' item addition, rerun."

"1" Item addition and rerun

If the sample is an unregistered sample, it is registered as a new request. If the sample is a registered item, the information of the registered item and item selection order from the host are subjected to the "OR" condition and then analyzed. If there is an order for an analyzed item, the item is analyzed again under rerun conditions. If there is an order for an uninspected item, the item is analyzed under initial-inspection conditions.

"2" No request

Applicable when there is no sample information on the host side or no new order for analysis. In this case, the request item box is omitted.

"3" Sample deletion

Sample registration is normally deleted from the system. Use this box when deleting orders by batch item registration from the host.

Item-selection instruction text – previous value exits

Applicable communication function: Batch-item registration (page 14)

Automatic item registration (page 16) Real item registration (page 18)

Communication direction: ADVIA 1650 Chemistry System ← Host Computer

| Maximum Frame Size | Maximum number of blocks |
|-----------------------|--------------------------|
| 256 bytes | 24 blocks |
| 512 bytes | 11 blocks |

Table 18 Item-selection Text - Previous Value Exists (first block, variable length)

| Data Item | Bytes | Code |
|---------------------------------|-------|--|
| Starting code | 1 | STX |
| Frame number | 1 | "1" – "7", "0" |
| Text classification | 1 | "o" (6FH), lowercase alphabetic oh |
| Equipment identification number | 1 | Not used, space (20H) |
| Total number of blocks | 2 | "01" to "24" (no zero suppression) |
| Block number | 2 | "01" (no zero suppression) |
| Number of items in a block | 3 | "001" to "999" (no zero suppression) |
| Sample classification | 1 | "N" (4EH) General sample "I" (49H) Interruption sample |
| Registration data | 1 | "0" (30H) New request. (If a registered sample [workorder] already exists, it is overwritten.) "1" (31H) Item addition, re-run. (In case of an unregistered sample, same as a new request.) "2" (32H) No request. (The ADVIA 1650 does not perform registration.) "3" (33H) Sample deletion. (Deletion of a previously registered sample) |
| Sample ID | 13 | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end Be sure to set unique sample numbers, even for an order of position No. designation. |

Table 18 Item-selection Text – Previous Value Exists (first block, variable length) (Cont)

| Data Item | Bytes | Code |
|--------------------------|---------------------|--|
| Position number | 7 | Filled with blank spaces when not used (20H). |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |
| Comment 1 | 16 | ASCII or shift JIS. Filled with blank spaces (20H) when not used. |
| Comment 2 | 16 | ASCII or shift JIS. Filled with blank spaces (20H) when not used. |
| Sex | 1 | "M" (4DH) male or "F" (46H) female: When not used, set to "M". |
| Age | 3 | 999 (right-justified): When not used, all three characters are 20H. Zero years old is 20H, 20H, "0". |
| Blood sampling date | 8 | YYYYMMDD Example: 19990229 Filled with blank spaces (20H) when not used. |
| Dilution coefficient | 4 | Format: 99.9 (right-justified): When not used, be sure to insert Format: ^1.0 (^0.1 - 99.9). ^ : blank space |
| Sample classification | 1 | "1" (31H) blood serum, "2" (32H) urine: When not used, be sure to input "1". |
| Container classification | 1 | "1" to "9" : When not used, "1". |
| Request item | (3+1+8+ | Item number 3 bytes, "999" (right-justified) |
| | 3) x n ¹ | Sample volume and dilution determined via the Analytical Parameters (Chemistry) window. One byte "M" or "D" or "U". Normally set to "M". |
| | | "M" = Analytical conditions values |
| | | "D" or "U" = Reanalysis condition values |
| | | Previous value: 8 bytes for 8-digit number (including "-" sign) and floating decimal point, right-justified ² Filled with blank spaces (20H) when not used. |
| | | Mark: 3 bytes Refer to mark specifications, page 36. |
| Spare | 1 | Not used (20H) |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF |

When the maximum frame size is 256 bytes, the maximum number of items is n = 11. When the maximum frame size is 512 bytes, the maximum number of items is n = 28.

$$123.45 \rightarrow ^{\land \land} 123.45$$

 $-6.7 \rightarrow ^{\land \land \land} -6.7$

² Previous values are as follows.

^{^:} Blank space (20H)

Table 19 Item-selection Text – Previous Value Exists (second to last block, variable length)

| Data Item | Bytes | Code |
|---------------------------------|--------------------|--|
| Starting code | 1 | STX |
| Frame number | 1 | "1" – "7", "0" |
| Text classification | 1 | "o" (6FH), lowercase alphabetic oh |
| Equipment identification number | 1 | Not used, space (20H) |
| Total number of blocks | 2 | "02" to "24" (no zero suppression) |
| Block number | 2 | "02" to "24" (no zero suppression) |
| Number of items in a block | 3 | "001" to "999" (no zero suppression) |
| Spare | 2 | Not used (20H) |
| Sample ID | 13 | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end |
| | | Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end |
| | | Be sure to set unique sample numbers, even for an order of position No. designation. |
| Position number | 7 | Filled with blank spaces when not used (20H). |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end |
| Request item | (3+1+8+ | Item number 3 bytes, "999" (right-justified) |
| | 3)x n ¹ | Sample volume and dilution determined via the Analytical Parameters (Chemistry) window. One byte "M" or "D" or "U". Normally set to "M". |
| | | "M" = Analytical conditions values |
| | | "D" or "U" = Reanalysis condition values |
| | | Previous value: 8 bytes for 8-digit number (including "-" sign) and floating decimal point, right-justified ² Filled with blank spaces (20H) when not used. |
| | | Mark: 3 bytes Refer to mark specifications, page 36. |
| Spare | 1 | Not used (20H) |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF |

When the maximum frame size is 256 bytes, the maximum number of items is n = 11. When the maximum frame size is 512 bytes, the maximum number of items is n = 28.

$$123.45 \rightarrow ^{\land \land} 123.45$$

 $-6.7 \rightarrow ^{\land \land \land} -6.7$

² Previous values are as follows.

^{^:} Blank space (20H)

Operational Notes:

- If the ID is classified as a sample ID for a batch-item inquiry or an item inquiry, set the sample ID in the itemselection instruction text, and leave the position number unused (barcode operation or creation of worksheets using the **ADVIA 1650** Chemistry System)
- If the ID is classified as a Sample Tray ID or rack-number ID for a batch-item inquiry or an item inquiry, set the position number to which the inquiry was addressed, and also set the sample number to be controlled by the host in the sample ID item (creation of worksheet using the host).
- At the start of Sample Tray analysis, the sample classified "I" as an interruption sample on sample classification is analyzed first. The ADVIA 1650 Chemistry System also has a built-in function to give priority to interruption samples for Sample Tray analysis. When the samples are analyzed on the rack option or conveyor-belt system, analysis begins sequentially as ordered regardless of that function. Therefore, interruption samples for these systems should be set on the Sample Tray.
- Registration data box

"0" New Request

A new request is unconditionally registered as a new item selection. In case of an already-registered sample, its information is deleted and overwritten, thereby erasing the analysis-result data. If you cannot judge whether it is a new request or registered item on the host side, specify all of the items in the next box as "1' item addition, rerun."

"1" Item addition and rerun

If the sample is an unregistered sample, it is registered as a new request. If the sample is a registered item, the information of the registered item and item selection order from the host are subjected to the "OR" condition and then analyzed. If there is an order for an analyzed item, the item is analyzed again under rerun conditions. If there is an order for an uninspected item, the item is analyzed under initial-inspection conditions.

"2" No request

Applicable when there is no sample information on the host side or no new order for analysis. In this case, the request item box is omitted.

"3" Sample deletion

Sample registration is normally deleted from the system. Use this box when deleting orders by batch item registration from the host.

Measurement Data Text

Batch data output (page 20) Real data output (page 21) Applicable communication function:

ADVIA 1650 Chemistry System → Host Computer Communication direction:

| Maximum Frame Size | Maximum number of blocks |
|--------------------|--------------------------|
| 256 bytes | 25 blocks |
| 512 bytes | 11 blocks |

Table 20 Measurement Data Text (first block, variable length)

| Data Item | Bytes | Code | | | |
|---------------------------------|-------|--|--|--|--|
| Starting code | 1 | STX | | | |
| Frame number | 1 | "1" – "7", "0" | | | |
| Text classification | 1 | "R" (52H) | | | |
| Equipment identification number | 1 | Not used, space (20H) | | | |
| Total number of blocks | 2 | "01" to "25" | | | |
| Block number | 2 | "01" (no zero suppression) | | | |
| Number of items in a block | 3 | "001" to "999" (no zero suppression) | | | |
| Inspection date | 8 | YYYYMMDD Example: 19990229 | | | |
| Sample classification | 1 | "N" (4EH) General sample "C" (43H) Control sample "S" (53H) STAT sample ¹ "I" (49H) Interruption sample ¹ | | | |
| ID specification | 1 | "0" (30H) New request. (If a registered sample [workorder] already exists, it is overwritten.) "1" (31H) Item addition, re-run. (In case of an unregistered sample, same as a new request.) "2" (32H) Rack number. (Not used.) "3" (33H) Sample deletion. (Deletion of a previously registered sample.) | | | |
| Sample ID | 13 | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end Be sure to set unique sample numbers, even for an order of potion No. designation. | | | |
| Position number | 7 | Filled with blank spaces when not used (20H). Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end | | | |

 Table 20 Measurement Data Text (first block, variable length) (Cont)

| Data Item | Bytes | Code | | | |
|--------------------------|--------------------------------|---|--|--|--|
| Comment 1 | 16 | ASCII or shift JIS. Filled with blank spaces (20H) when not use | | | |
| Comment 2 | 16 | ASCII or shift JIS. Filled with blank spaces (20H) when not used. | | | |
| Sex | 1 | "M" (4DH) male or "F" (46H) female: When not used, set it to "M". | | | |
| Age | 3 | 999 (right-justified): When not used, age is always "0." | | | |
| Blood sampling date | 8 | YYYYMMDD Example: 19990229 | | | |
| Dilution coefficient | 4 | Format: 99.9 (right-justified): When not used, be sure to insert Format: ^1.0 (^0.1 - 99.9). ^ : blank space | | | |
| Sample classification | 1 | "1" (31H) blood serum, "2" (32H) urine: When not used, be sure to set "1". | | | |
| Container classification | 1 | "1" to "9" : When not used, "1". | | | |
| Request item | (3+1+8+ 3) x n ² | Item number 3 bytes, "999" (right-justified) | | | |
| | | Analysis condition 1 byte "M" or "D" or "U" Normally set to "M" | | | |
| | | Previous value: 8 bytes for 8-digit number (including "-" sign) and floating decimal point, right-justified Filled with blank spaces (20H) when not used. | | | |
| | | Mark: 3 bytes Refer to mark specifications, page 36. | | | |
| Spare | 1 | Not used (20H) | | | |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF | | | |

STAT samples and control samples are sent with the sample classification and a unique alphanumeric ID assigned by the ADVIA 1650.

Table 21 Measurement Data text (second to last block, variable length)

| Data Item | Bytes | Code | | |
|---------------------------------|-------|---|--|--|
| Starting code | 1 | STX | | |
| Frame number | 1 | "1" – "7", "0" | | |
| Text classification | 1 | "R" (52H) | | |
| Equipment identification number | 1 | Not used, space (20H) | | |
| Total number of blocks | 2 | "02" to "25" (no zero suppression) | | |
| Block number | 2 | "02" to "25" (no zero suppression) | | |
| Number of items in a block | 3 | "001" to "999" (no zero suppression) | | |
| Inspection date | 8 | YYYYMMDD Example: 19990229 | | |
| Sample classification | 1 | "N" (4EH) General sample "C" (43H) Control sample "S" (53H) STAT sample "I" (49H) Interruption sample | | |

STAT ID: 3 characters (E01 to E84), no zero suppression, left justified, 20H at end.

Control ID: 4 characters (PA01 to PJ 20) no zero suppression, left justified, 20H at end.

When the maximum frame size is 256 bytes, the maximum number of items is n = 10. When the maximum frame size is 512 bytes, the maximum number of items is n = 27.

Table 21 Measurement Data text (second to last block, variable length) (Cont)

| Data Item | Bytes | Code | | | | |
|------------------|-------------------------------|--|--|--|--|--|
| ID specification | 1 | "0" (30H) New request. (If a registered sample [workorder] already exists, it is overwritten.) | | | | |
| | | "1" (31H) Item addition, re-run. (In case of an unregistered sample, same as a new request.) | | | | |
| | | "2" (32H) Change attribute. (The system does not perform registration.) | | | | |
| | | "3" (33H) Sample deletion. (Deletion of a previously registered sample.) | | | | |
| Sample ID | 13 | Arbitrary-sample ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end | | | | |
| | | Barcode ID: Up to 13 alphanumeric characters, no zero suppression, left-justified, 20H at the end | | | | |
| | | Sequential ID: Up to 9 digits, no zero suppression, left-justified, 20H at the end | | | | |
| | | Be sure to set unique sample numbers, even for an order of position No. designation. | | | | |
| Position number | 7 | Filled with blank spaces when not used (20H). | | | | |
| | | Sample Tray: Fixed-length, 5 digits, no zero suppression, left-justified, 20H at the end | | | | |
| | | Rack No.: Fixed-length, 7 digits, no zero suppression, left-justified | | | | |
| | | Independent holder No.: Fixed-length, 4 digits, no zero suppression, left-justified, 20H at the end | | | | |
| Request item | (3+1+8+3) x n ¹ | Item number 3 bytes, "999" (right-justified) | | | | |
| | | Analysis condition 1 byte "M" or "D" or "U". Normally set to "M". | | | | |
| | | Previous value: 8 bytes for 8-digit number (including "-" sign) and floating decimal point, right-justified ² Filled with blank spaces (20H) when not used. | | | | |
| | | Mark: 3 bytes Refer to mark specifications, page 36. | | | | |
| Spare | 1 | Not used (20H) | | | | |
| Termination code | 5 | (ETX or ETB) + checksum + CR + LF | | | | |

When the maximum frame size is 256 bytes, the maximum number of items is n = 13. When the maximum frame size is 512 bytes, the maximum number of items is n = 31.

$$\begin{array}{c} 123.45 \rightarrow \text{^{^{^{^{^{}}}}}}123.45 \\ -6.7 \rightarrow \text{^{^{^{^{^{^{}}}}}}}-6.7 \\ \text{Overflow} \rightarrow \text{////////} \end{array}$$

No measurement value → ^^^^^^

^: Blank space (20H)

Items that are set to output qualitative measurement data are as follows:

- Each of the items is fixed at 8 characters in length. Up to 6 characters are used for analysis-condition settings and they are right-justified. The first 2 characters are used for spaces. Setting 3 two-byte characters is possible, and in this case, SHIFT JIS codes will be output.
- The position-number box of the sample that was subjected to barcode analysis using a sample ID is invalid.

Result formats are as follows.

Mark (Result Flag) Specifications
A mark or result flag consists of a total of 3 bytes: judgment, status, and rerun, each 1 byte.

| Judgment (first byte) | | | Status (seco | Down (third by to) | | | |
|-----------------------|--------------------------|-------|-------------------------------------|--------------------|--|--------------------|-----------------|
| | | Pl | notometric Flag | Electrolyte Flag | | Rerun (third byte) | |
| С | Calibration not possible | s | Insufficient sample | S | Insufficient sample | R | Rerun value |
| L, H | Abnormal value limit | t | Insufficient diluent liquid | S | Safety error | Space | first run value |
| I, h | Normal value limit | S | Safety | В | Liquid remain- ing in dilution cup | | |
| Space | No mark (flag) | r | Insufficient reagent | Т | Abnormal thermistor | | |
| | | u, d | Light absorption limit | r | Insufficient reagent | | |
| | | р | Prozone limit | u | Selectivity abnormal | | |
| | | * | Dispersion abnor- mality | * | CAL allow- able value | | |
| | | n | Abnormal number of effective points | Space | No mark (flag) | | |
| | | N | Cell-blank abnormality | | | | |
| | | U, D | Light absorption limit | | | | |
| | | Space | Space No mark (flag) | | | | |

In the above table, the higher the position in the mark (result flag) column, the higher the degree of priority.