



The List of Test Result Transmission Regarding to The Sample Kit

(Accompany With Bioneer ExiStation™ LIS Interface Specification)

Document Number: BIS-JG-104-1-1.0.0.1

List of Test Result Transmission

The Test Result of each sample has IPC CT, IPC Result, CT, Delta Rn, Result and etc. ExiStation™ transmits such results above; These result items are transmitted individually in Test Result Records.

Due to the item that will be transmitted is dependent on a Sample Kit(s), please refer to the list of result items in the table provided below.

* The following notation, [Sample Kit Name], is replaced with the name of Bioneer Sample kit.

* The following notation of [1st Target Name] or [2nd Target Name] is replaced with the name of the Target.

Classification	Sample Kit Name	Sending Result Items
Quality	HAV-QL-ES	1. IPC CT
	IFB-QL-ES	2. IPC Result
	MTB-QL-ES	3. [Sample Kit Name] CT
	CMV-QL-ES	4. [Sample Kit Name] Delta Rn
	CPN-QL-ES	5. [Sample Kit Name] Result
	EBV-QL-ES	
	MPN-QL-ES	
2-Quality	CNM-QL-ES	1. IPC CT
		2. IPC Result
	SAM-QL-ES	3. [1st Target Name] CT
		4. [1st Target Name] Delta Rn
2-Quality	HPV-QL-ES	5. [1st Target Name] Result
		6. [2nd Target Name] CT
	MTN-QL-ES	7. [2nd Target Name] Delta Rn
		8. [2nd Target Name] Result
Quantity	CMV-QN-ES	1. IPC CT
	EBV-QN-ES	2. IPC Result

	HBV-QN-ES	3. [Sample Kit Name] CT
	HCV-QN-ES	4. [Sample Kit Name] (copy/rxn)
	HIV1-QN-ES	5. [Sample Kit Name] (copy/ml)
	HIV2-QN-ES	6. [Sample Kit Name] (IU/ml) 7. [Sample Kit Name] Result
2-Quantity	HIM-QN-ES	1. IPC CT 2. IPC Result 3. [1st Target Name] CT 4. [1st Target Name] (copy/rxn) 5. [1st Target Name] (copy/ml) 6. [1st Target Name] (IU/ml) 7. [1st Target Name] Result 8. [2nd Target Name] CT 9. [2nd Target Name] (copy/rxn) 10. [2nd Target Name] (copy/ml) 11. [2nd Target Name] (IU/ml) 12. [2nd Target Name] Result